


OASYS®

Surgical Technique

Occipito-Cervico-Thoracic System



System Overview

The OASYS® Occipito-Cervico-Thoracic System was developed to provide the surgeon with versatility for the treatment of pathologies of the occipitocervical junction, and the posterior cervical and upper thoracic spine. The modular system consists of hooks, polyaxial screws and rods as well as plates and bone screws. A variety of occiput plates and connectors complement the product line.





Table of Contents

System Overview	2
Implant Overview	4
Surgical Procedure	10
Polyaxial Screw Placement	10
Rod Placement	13
Blocker Insertion	15
Transverse Connection	17
Offset & Rod to Rod Connectors	19
Hook Placement	20
Occipital Fixation	21
Midline Occiput Plate Placement	23
Bilateral Occiput Plate Placement	30
Implants	33
Instruments	36



Stryker® Spine would like to thank
John J. Carbone, MD
Baltimore, Maryland
for his contribution.

Implant Overview

Polyaxial Screws

One of the key features of the Biased Angle Polyaxial Screw is the offset angle of the screw head, which allows for combined **divergent screw angulation of 110° (or 55° in one direction)**.

The high degree of angulation promotes screw placement at an optimal anatomic position, helping to simplify the surgical procedure by minimizing the need for rod contouring due to linear placement of the screw heads.


There are five types of polyaxial screws available in the OASYS[®] System. (See complete screw listing on page 33-34).

The standard **Biased Angle Screws** are available in diameters of 3.5mm (10 to 54mm lengths in 2mm increments) and 4.0mm (10 to 42mm in 2mm increments, 46 and 50mm lengths). The 4.0mm Biased Angle Screws are anodized blue for easy identification.


The **Medial Biased Angle Screws** are available in diameters of 3.5mm (20, 22, and 24 to 40mm lengths in 4mm increments) and 4.0mm diameters (20 to 52mm lengths in 4mm increments). The tulip heads of the Medial Biased Angle Screws are anodized green for easy identification.

The **Smooth Shank Biased Angle Screws** are available in a 3.5mm diameter (22 to 40mm lengths in 2mm increments). These screws have a 10mm non-threaded shank to potentially avoid nerve root or tissue irritation. The tulip heads of the Smooth Shank Biased Angle Screws are anodized blue for easy identification.

The **Cancellous Biased Angle Screws** are available in diameters of 3.5mm and 4.0mm (10 to 24mm lengths in 2mm increments). These screws have a thread geometry designed for purchase in cancellous bone. The tulip heads of the Cancellous Biased Angle Screws are anodized blue-green for easy identification.



Medial Biased Angle Screw



The **Non-Biased Angle Polyaxial Screws** are available in diameters of 3.5mm, 4.0mm and 4.5mm. The 3.5mm and 4.0mm screws are 10 to 20mm in length in 2mm increments while the 4.5mm screws are 20 to 52mm in length in 4mm increments. The Non-Biased Angle Polyaxial Screws offer up to 30° of angulation in each direction or 60° conical. These screws can be used in instances where the surgeon does not need the angulation offered with the Biased Angle Screws. A combined color approach has been used to identify the screw type and diameter.



Non-Biased Angle Polyaxial Screw

The tulip heads of the Non-Biased Angle Polyaxial Screws are anodized fuchsia while the 4.5mm bone screws are anodized purple.

Screw Type	Diameter	Lengths	Tulip Head Color	
Biased Angle Screws	3.5mm	10 - 54mm (in 2mm increments)	Silver	
	4.0mm	10 - 42mm, 46 and 50mm (in 2mm increments)	Silver	
Medial Biased Angle Screws	3.5mm	20, 22 and 24 - 40mm (in 4mm increments)	Green	
	4.0mm	20 - 52mm (in 4mm increments)	Green	
Smooth Shank Biased Angle Screws	3.5mm	22 - 40mm (in 2mm increments)	Blue	
Cancellous Biased Angle Screws	3.5mm	10 - 24mm (in 2mm increments)	Blue-Green	
	4.0mm	10 - 24mm (in 2mm increments)	Blue-Green	
Non-Biased Angle Polyaxial Screws	3.5mm	10 - 20mm (in 2mm increments)	Fuchsia	
	4.0mm	10 - 20mm (in 2mm increments)	Fuchsia	
	4.5mm	20 - 52mm (in 4mm increments)	Fuchsia	

Hooks

A variety of laminar **Hooks** are available to better match the individual patient's anatomy. Standard hooks are available in four throat heights: 3.5, 5.0, 6.5 and 8.0mm. Right and left offset hooks are also available for both right and left application.



Part #	Description
48551049	Hook, Standard, 3.5mm
48551050	Hook, Standard Short, 5.0mm
48551053	Hook, Standard, 6.5mm
48551055	Hook, Standard Tall, 8.0mm
48551060	Hook, Offset Right
48551065	Hook, Offset Left



Hook, Standard, 6.5mm



Hook, Standard Tall, 8.0mm



Hook, Standard, 3.5mm



Hook, Standard Short, 5.0mm



Hook, Offset Right



Hook, Offset Left

Rods

A variety of **Rod** offerings are available to provide intraoperative solutions. For added stiffness, Vitallium[®] is available. Titanium and Vitallium[®] pre-cut offerings are available with the OASYS System, as well as Titanium Commercially Pure in 240mm.

Part #	Description
48553080	3.5mm x 80mm Vitallium [®] Rod
48553120	3.5mm x 120mm Vitallium [®] Rod
48553240	3.5mm x 240mm Vitallium [®] Rod
48553350	3.5mm x 350mm Vitallium [®] Rod

48552025	3.5mm x 25mm Titanium Alloy Rod
48552030	3.5mm x 30mm Titanium Alloy Rod
48552040	3.5mm x 40mm Titanium Alloy Rod
48552050	3.5mm x 50mm Titanium Alloy Rod
48552060	3.5mm x 60mm Titanium Alloy Rod
48552070	3.5mm x 70mm Titanium Alloy Rod
48552080	3.5mm x 80mm Titanium Alloy Rod
48552120	3.5mm x 120mm Titanium Alloy Rod
48552240	3.5mm x 240mm Titanium Alloy Rod
48551240	3.5mm x 240mm Titanium, Commercially Pure Rod
48551350	3.5mm x 350mm Titanium Alloy Rod



Occipital Fixation




The OASYS® System offers a midline occiput plate option as well as bilateral occiput plates to better match individual patient anatomy.

Part #	Description
48551040	130° Occiput Plate
48551041L	130° Pre-contoured Occiput Plate, Left
48551041R	130° Pre-contoured Occiput Plate, Right
48551042	100° Occiput Plate
48551043L	100° Pre-contoured Occiput Plate, Left
48551043R	100° Pre-contoured Occiput Plate, Right
48551044	Midline Occiput Plate, Small
48551045	Midline Occiput Plate, Medium
48551046	Midline Occiput Plate, Large
48551047	Midline Occiput Plate, Large Long
48551048	Midline Occiput Plate, Mini



Bone Screws

The following bone screws are available for use with the OASYS® plates:

Diameter	Lengths (in 2mm increments)	Color	
3.5mm	6 – 16mm	Silver	
4.0mm	6 – 16mm	Blue	
4.5mm	6 – 16mm	Purple	


Connectors

Several types of **Connectors** are available in the OASYS[®] System.

Transverse Connector

The **Transverse Connector** is designed for added torsional rigidity of a bilateral rod construct. The clip-on design features a very low profile and does not add significantly to the height of the construct. The transverse bar is available in 40, 60, and 80mm lengths. The bar of the transverse connector can be cut so that the length can be adjusted to better match individual patient anatomy.

Part #	Description
48551070	Transverse Connector, 80mm
48551071	Transverse Connector, 60mm
48551072	Transverse Connector, 40mm



Transverse Connector

Cross Connector Plate

In cases where screws are closer together, this arched implant allows connection to the tulip heads. The arch is designed for additional clearance above the dura and spinal cord. The **Cross Connector Plate** is available in 24, 32, and 40mm lengths with 10mm of additional variability for each plate. The plates can be bent to adapt to variable patient anatomy.

Part #	Description
48551073	Cross Connector Plate, 24mm
48551074	Cross Connector Plate, 32mm
48551075	Cross Connector Plate, 40mm



Cross Connector Plate

Offset Connector

As an alternative to offset rod contouring in cases of non-linear screw or hook placement, **Offset Connectors** can be used to provide additional medial-lateral offset. The bar is available in 12 and 20mm lengths.

Part #	Description
48551080	Offset Connector, 20mm
48551081	Offset Connector, 12 mm



Offset Connector

Rod-to-Rod Connector

The **Rod-to-Rod Connector** facilitates connection of the 3.5mm diameter rod to a spine construct using:

- A 6.0mm diameter rod: Xia®
- A 5.5mm diameter SR90D®
- A 4.5mm diameter Xia® rod
- Another 3.5mm diameter OASYS® rod

Part #	Description
48551088	3.5mm to 3.5mm Parallel Rod-to-Rod Connector
48551091	3.5mm to 4.5mm Parallel Rod-to-Rod Connector
48551089	3.5mm to 5.5mm Parallel Rod-to-Rod Connector
48551090	3.5mm to 6.0mm Parallel Rod-to-Rod Connector
48551085	3.5mm to 3.5mm Axial Connector
48551084	3.5mm to 4.5mm Axial Connector
48551086	3.5mm to 5.5mm Axial Connector
48551087	3.5mm to 6.0mm Axial Connector




Parallel Rod-to-Rod Connector



Axial Connector

Saddle Connector

The **Saddle Connector** is an additional connector option which helps to allow connection of the 3.5mm rod to a spine construct using Xia® II and/or Xia®3* Polyaxial Screws. The Saddle Connector is to allow the use of a larger size bone screw in the thoracic area of T1 to T3 per the current OASYS® indications.



Saddle Connector
48551094

*Refer to the Xia® Surgical Technique Guide for information on screw insertion and blocker final tightening.

Surgical Procedure

Preoperative Planning

Use the appropriate imaging techniques to outline the patient's osseous anatomy and to determine the proper size and type of instrumentation to be used. Identify the components to be used for the assembly. Keep in mind that changes to the final configuration may become necessary based on intraoperative findings.

Patient Positioning and Exposure

The patient is placed in the prone position with head and neck held securely in optimum alignment. A standard midline incision is performed at the appropriate levels. The exposure may be extended for one or two levels below the inferior end of the planned fusion to allow for easy placement of the instrumentation.

Polyaxial Screw Placement

Following exposure, penetrate the cortex with an **Awl**, burr, or a drill to mark the entry point of all the screws using anatomic, fluoroscopic, or image-guided technique. To minimize the need for rod contouring and for easy rod insertion, it is advisable to align the screw holes as much as possible.

In applicable anatomic areas, a **Pedicle Probe** – followed by a **Tap** – can be used to prepare the screw pathway. In all other cases, the preparation of the screw hole should follow the steps below, using the appropriate Drill Bit and Tap.

The **Adjustable Drill Guide** allows for a single **Drill Bit** to be used for preparation of variable depths. To set the depth which corresponds to the final screw length, grasp knob, then rotate the locking pin into the slot corresponding to the desired depth. Push the pin forward and lock the fixation nut at the end of the sleeve by rotating clockwise.



Awl
48560010



Adjustable Drill Guide
48561011

If the sleeve is not moving freely, loosen the fixation nut before manipulating the pin. Drill depth can be adjusted while the Drill Bit is inserted in the Adjustable Drill Guide, however, it is not recommended while drilling. The fixation nut and the inner sleeve may be disassembled for cleaning purposes.

For common screw lengths, **Fixed Drill Guides** are also available. These Fixed Drill Guides are available for 12 and 14mm screws. The 14mm guide has a gold band around it for easy identification.



Fixed Drill Guide

Select the **Drill Bit** that is appropriate for the size of the screw to be used (3.5mm, 4.0mm, or 4.5mm). The Drill Bits are undersized by 1mm (i.e., the 3.5mm drill has a diameter of 2.5mm), corresponding to the inner diameter of the screws. The Drill Bits can be attached to the **Quick-Release Handle** or used with standard power equipment (AO attachment).

Part #	Description
48561034	Fixed Drill Guide, 12mm
48561035	Fixed Drill Guide, 14mm

Part #	Description
48560323	2.5mm Drill for 3.5mm Polyaxial Screw
48560423	3.0mm Drill for 4.0 & 4.5mm Polyaxial Screws

It is recommended that a Tap be used to finalize the preparation of the screw pathway. Taps are available in three sizes, 3.5mm, 4.0mm, and 4.5mm. There are two types of taps - cancellous and cortical - to match the thread pitch of the screws being used. The optional **Tap Sleeve** may be used to accurately measure the tapping depth or to protect soft tissue.



Tap Sleeve
48561914
48561915

Note: Using incorrect Tap or Drill combination can result in reduced bone purchase or possibility of instrument breakage.



Cortical Tap

Part #	Description
48560314	3.5mm Tap
48560414	4.0mm Tap
48560514	4.5mm Tap



Cancellous Tap

Part #	Description
48561916	3.5mm Cancellous Tap
48561917	4.0mm Cancellous Tap

The **Ball Tip Probe** can be used to feel the integrity of the pedicle wall after drilling or tapping. The probe has a 1.8mm ball tip for smaller pedicles. The **Depth Gauge** can be used to confirm the screw length.

Note: The 4.0mm Drill followed by the 4.5mm Tap should be used to prepare the pathway for the 4.5mm Non-Biased Angle Screws. Also, when preparing the pathway for the Cancellous Screws, make sure to use the Cancellous Tap. The Cancellous Tap has a different thread geometry to match the geometry of the Cancellous Screws. The standard Tap Sleeve can be used with both the 3.5mm and 4.0mm Standard and Cancellous Taps. For the 4.5mm Tap, use the 4.5mm Tap Sleeve.

A polyaxial screw is attached to the **Polyaxial Screwdriver** by placing the hex head socket over the bone screw and threading the end of the outer sleeve into the polyaxial screw head for stable fixation. Prior to screw implantation, use the gauges incorporated into the implant tray to verify the screw dimensions. Aligning the Polyaxial Screwdriver in the same axis as the screw hole facilitates screw insertion. Once the screw has been inserted, the Polyaxial Screwdriver can be removed by unthreading the end of the sleeve from the screw. There is an optional outer sleeve that may be used on the Polyaxial Screwdriver. This sleeve allows you to hold the shaft of the Polyaxial Screwdriver for better stabilization without having the polyaxial screw disengage from the Polyaxial Screwdriver. The sleeve can be placed on the Polyaxial Screwdriver by holding in the button while sliding the sleeve over the shaft.

Tip: If the tip of a biased angle screw hits the bone before the screw is fully inserted, and the screw is unable to be driven in further, unthread the screwdriver sleeve from the tulip head while leaving the hex socket engaged with the screw to allow the screw head to spin freely and the bone screw to be driven in further.



Tap with Quick Release Handle
48560012



Tap Sleeve
48561914
48561915



Polyaxial Screwdriver
48562016



Note: Ensure the screwdriver sleeve is tightly assembled to the Polyaxial Screwdriver shaft before attempting to attach a screw.

Note: To clean the Polyaxial Screwdriver, remove the sleeve by pressing the button. The holes in the shaft of the Polyaxial Screwdriver are designed to allow flushing of the instrument.


To minimize the potential for rigid locking and limiting the range of angulation, do not seat the screw head too tightly to the bone. Should the screw head become locked to the bone screw, use the **Screw Head Adjuster** to break the lock. If needed, the screw may be backed out about half a turn to allow for the head to move freely.

If it becomes necessary to remove a screw, begin by unthreading the bone screw with the hex portion of the Polyaxial Screwdriver, and then reattach the outer sleeve to complete the removal.

For additional height adjustments to the screw, the **Poly Adjustment Driver** can be used in place of the Polyaxial Screwdriver. This instrument has the same internal hex as the Polyaxial Screwdriver, but has no sleeve to impede vision.

Rod Placement

Once all the bone screws have been inserted and aligned, a **Rod Template** is available to estimate the required rod length and contour. If the screw heads need to be aligned in a more linear fashion, the Screw Head Adjuster can be used for this purpose. Because of the angle of the biased cut and the corresponding range of angulation of the screw, it is recommended that the notched lines on the screw head face in the direction in which the maximum bone screw angulation is desired (corresponds to the smallest angle between the screw and the rod). All screws should be aligned prior to rod and blocker insertion as further adjustments will become difficult once the construct has been assembled.



**Screw Head Adjuster
48560021**



**Rod Template
48560017**

To prepare the **Rod/Plate Cutter** for rod insertion, move the handles into the “open” position as indicated with the laser-marked arrows. Insert the rod into the head of the Rod/Plate Cutter from the side with the laser-marked arrows, and compress the handles as indicated. Since the rod will be cut slightly below the surface of the cutting head, there will be approximately 3mm added to the rod length from the point of insertion. It is recommended that the rod be cut to the desired length prior to contouring as it may become difficult to cut a rod that has already been bent.

Note: Lubrication of the Rod/Plate Cutter is required to maintain ease of use of the instrument.

A simple, “pin cutter” style **Small Rod Cutter** is also available to cut the rods. This cutter is for 3.5mm diameter rods only.



Rod/Plate Cutter
48560018



Small Rod Cutter
48561038

Using the **Rod/Plate Bender** (the side with the round bending iron), the rod can be contoured to match the sagittal alignment of the spine as well as the coronal orientation of the screws. Avoid sharp, excessive, or repeated bending of the rod to maintain material integrity. Prior to implantation, inspect the rod for any damage (e.g., notching) it may have sustained during preparation.

If damaged, use new rod. **In-situ Benders** are also available for additional contour adjustments.

Once the rod has been shaped to its final configuration, introduce the rod into the screw heads using the **Rod Forceps**.

The **Rod Rotation Forceps** allow for minor rod rotation maneuvers. The tips have serrated teeth for better grip on the rod.



Rod Plate Bender
48560019



In-situ Benders
48560022 - Right
48560023 - Left



Rod Rotation Forceps
48561039

Blocker Insertion

There are four options for insertion and alignment of the Blockers:

1. Locking Persuader and Hexdriver
2. Persuader and Hexdriver
3. Insertion Tube and Hexdriver
4. Hexdriver by itself


To facilitate rod insertion, a **Persuader** may be applied to fully seat the rod into the screw head. While keeping the instrument aligned with the course of the rod, slide the tip of the Persuader over the screw head. By gently squeezing the trigger, the rod will be depressed into the screw head. The shaft of the Persuader is cannulated to allow for the passage of the **Hexdriver** and engagement of the blocker while the Persuader is in place.



**Persuader
48560024**

Note: The Persuader is not used to cause spinal correction in patient.

The **Locking Persuader** has two bars that push down on both sides of the rod. It also has a ratcheted trigger design that locks the Locking Persuader arms and allows you to maintain the rod reduction without requiring a tight grip on the trigger. When attaching to the screw head, align the side windows parallel to the rod. To disengage, press the thumb release button, rotate the Locking Persuader slightly to either side, and gently pull it off the screw head.



**Locking Persuader
48561024**

When using the non-locking Persuader, you must maintain pressure on the trigger while the blocker is inserted. To disengage, release the trigger, rotate the Persuader slightly to either side, and gently pull it off the screw head.

As an alternative to using the Persuader, the cannulated **Insertion Tube/Anti-Torque Key** can be used to push the rod down into the screw head and align the Polyaxial Screwdriver during blocker insertion. A specially designed ring on the inside of the **Insertion Tube** provides an interference fit with a fully inserted Hexdriver to potentially reduce the risk for an accidental disassembly. It is recommended to hold both instruments during their removal from the rod construct.




**Insertion Tube/Anti-Torque Key
48560027**

Load the blocker firmly onto the tip of the **Hexdriver** and place it into the screw head, using the cannulated instruments as guides or in a free-hand manner. For proper loading onto the Hexdriver, the blockers should be stored in the tray with the laser-marked line facing up. The buttress thread design of the blocker and proper alignment of the Hexdriver minimize the potential for cross-threading.* Turning the blocker counter-clockwise about a quarter turn when first introduced facilitates the insertion into the screw head or hook. When used with the Persuader or with the Insertion Tube, laser markings on the Hexdriver aid in visually assessing the depth of blocker insertion. The lower of the two lines indicates the level at which the blocker engages the screw head; the Hexdriver should be removed when the upper laser marking has been reached.



Note: The Hexdriver should be used for blocker insertion and provisional tightening only. Do not perform final tightening with the Hexdriver as it will result in damage to the instrument over time and incorrect final implant torque.

Once all of the blockers have been inserted, the Insertion Tube or the Persuader, in combination with the **Torque Wrench**, should be used for final tightening. The blocker is completely tightened when the two arrows on the shaft of the Torque Wrench are aligned, which corresponds to 3Nm.



Note: It is important to tighten the blocker to the recommended torque to ensure future integrity of the construct. Tightening under or over the torque limit is inadvisable and should be avoided. Appropriate counter-torque must be applied with the Insertion Tube.



Torque Wrench
48560028



Insertion Tube/Anti-Torque Key
48560027

Should compression of the adjacent instrumented levels be desired prior to final blocker tightening, it can be achieved by placing the **Compressor** over the neighboring screw heads or hooks and using the ratchet mechanism to adjust the amount of compression needed. Similarly, the **Distractor** may be used to achieve distraction of adjacent levels.



Compressor
48561026



Distractor
48561041

*Data on file at Stryker Spine

Transverse Connection

Two options are available for transverse connection to offer additional torsional rigidity for bilateral constructs:

1. Transverse Connectors
2. Cross Connector Plates

The **Transverse Connector** attaches directly to the rod resulting in a low profile construct. To insert, measure the distance between the rods, cut the connector bar with the Rod/Plate Cutter to the appropriate length (distance between the rods plus 10mm for adequate fit within the connector clip), and pre-load both connector clips onto the connector bar. If the bar is contoured with the Rod/Plate Bender, use the appropriate precautions as described previously.


Note: Usage of the Connector Clip Inserter is strongly recommended to neutralize the forces needed for assembly and prevent damage to implant.

Place the hooked part of the Connector Clip Inserter under the rod, position the connector on the pin protruding from the sliding upper shaft, and snap the clip onto the rod by squeezing the Connector Clip Inserter handle. Repeat the insertion steps for the contra-lateral clip. Check that the connector bar protrudes outside both clips for proper locking. To secure both clips in their final position, the Insertion Tube in combination with the Torque Wrench should be used for final tightening. The recommended torque value is the same as for the blockers: 3Nm.

If the Transverse Connector needs to be removed from the rod construct, remove the set screws from both clips using the Hexdriver and pull out the transverse bar. By advancing the **Connector Clip Remover** into the head of the clip, the clip will be released from the rod.



Transverse Connector



Connector Clip Inserter
48560029



Connector Clip Inserter



Connector Clip Remover
48560129

The **Cross Connector Plate** attaches to the tulip heads to accommodate narrow spaces between the screws. When using this type of connector, you must use a connector plate blocker instead of the standard blocker to tighten the rod to the two screws where you wish to make the transverse connection. A locking nut is then used to tighten the plate onto the identified tulip heads.


After deciding where to add the cross connection, insert the standard blockers on all screws except where the connection is desired. Insert the connector plate blockers onto those tulip heads using the Hexdriver.

Note: The Persuaders do not fit over the connector plate blockers. Thus, if persuasion is needed, the Persuader or Insertion Tube must be attached to the tulip heads of the levels adjacent to the cross connector placement.


Once the rod placement is verified, perform final tightening on all blockers. For final tightening of the connector plate blockers, place the Torque Wrench on the blocker and place the **Anti-Torque Key** on an adjacent level.

Next, choose the appropriate size Connector Plate (24, 32, or 40mm) by measuring the distance between the rods. The Connector Plate's screw hole is slotted to accommodate up to an additional 10mm of distance between the rods. You can further contour the Connector Plate using the Connector Plate Bender. To bend the Connector Plate, slide one side of the plate fully into the first bender (benders are exactly the same—there is no left or right). Then slide the other side of the plate into the second bender. To straighten the arch, insert the plate with the arch facing down. To increase the bend, insert the plate with the arch facing up. You can also apply a slight twist as needed.

Note: To maximize strength of the construct, it is important that the Connector Plate lay flatly on the blockers. If the Connector Plate does not fit properly, contour the plate using the Connector Plate Benders.



Extended blocker placed onto medial-biased tulip head.




Anti-Torque Key is placed at adjacent level. Final tighten extended blockers to 3Nm with Torque Wrench.



Bend Cross Connector Plate, if necessary.

Next, lay the Connector Plate on the connector plate blockers. Attach the **Locking Nut Socket** to the Torque Wrench, and use the assembled instrument to place the locking nut over the Connector Plate. After verifying the correct fit of the Connector Plate, tighten the locking nuts by attaching the Anti-Torque Key to the adjacent tulip head and aligning the Torque Wrench arrows to achieve the 3Nm torque.




Attach Locking Nut Socket to Torque Wrench.

Offset and Rod-to-Rod Connectors


In cases in which the patient's anatomy requires significantly different lateral or medial screw positions, Offset Connectors may be used to facilitate rod attachment.

The bar of the **Offset Connector** is available in two sizes: 12mm and 20mm. The 12mm length allows for up to 4mm of offset from the axial alignment, while the 20mm connector will allow for up to 12mm. If desired, the bar can be shortened by using the Rod/Plate Cutter as described previously.



Offset Connector

The **Rod-to-Rod Connectors** are available to connect the system to a construct with a 6.0, 5.5, 4.5, or 3.5mm diameter rod. All four sizes are available in a side-to-side and axial version. It is recommended that the Rod-to-Rod Connector be preloaded on the larger rod: Xia[®] or SR90D[®]. Slide the 3.5mm rod into the connector, and then use the Hexdriver to provisionally tighten the connector set screws to secure the assembly in place.



Axial Connector

The Torque Wrench must be used for final tightening of all connector set screws. The recommended torque value is the same as for the blockers: 3Nm.

Note: The connector set screws used in the Transverse, Offset, and Rod-to-Rod Connectors are smaller than the blockers used in combination with the Polyaxial Screws, Hooks, and the Midline Occiput Plates, and they are not interchangeable.




Parallel Rod-to-Rod Connector

Saddle Connector

The **Saddle Connector** is available to connect the OASYS[®] Occipital-Cervico-Thoracic System to a construct with Xia[®] II or Xia[®]3 Polyaxial screws.* The Saddle Connector is to allow the use of a larger size bone screw in the thoracic area of T1 to T3 per the current OASYS[®] indications.

To use the Saddle Connector, first follow the surgical technique for screw insertion for the Xia[®] II or Xia[®] 3 polyaxial screws and the Oasys[®] surgical technique for placement of Oasys[®] polyaxial screws. Place the Oasys[®] Saddle Connector in the tulip-head of the Xia[®] screw. The Saddle Connector should be facing upwards in a U-shape. Introduce the Oasys[®] rod into the Saddle Connector of the Xia[®] screw and into the adjacent Oasys[®] tulip heads. Fix the rod to the screw using standard Xia[®] blockers in the Xia[®] screw tulip heads and standard Oasys[®] blockers in the Oasys[®] screw tulip heads.




Saddle Connector in Polyaxial Screw

Hook Placement

Lamina preparation and hook size templating is facilitated by using the **Hook Preparer**. The two blades of the Hook Preparer – one on each end – correspond to the short and tall standard hook blades.

Select the appropriate Hook and use the **Hook Forceps** to place it into position.

For rod sizing, contouring, and insertion as well as final construct fixation (i.e., blocker insertion), follow the steps described in the section on Polyaxial Screws above.



**Hook Preparer
48560032**



**Hook Forceps
48560033**


*Refer to the Xia[®] Surgical Technique Guide for information on screw insertion and blocker final tightening.

Occipital Fixation

Implant Overview

Bone Screws


The following bone screws are available for use with the OASYS® plates:

Diameter	Lengths (in 2mm increments)	Color	
3.5mm	6 – 16mm	Silver	
4.0mm	6 – 16mm	Blue	
4.5mm	6 – 16mm	Purple	

Midline Occiput Plates

The OASYS® system features a Midline Occiput Plate. This single plate allows for positioning directly on the midline keel of the occiput.

The Midline Occiput Plate is uniquely designed to contour to the patient's anatomy to reduce the need or amount of rod-contouring. The plate is available in four standard sizes*: Small, Medium, Large and Large Long. It can be described as having a plate **body** that attaches to the occiput, two **legs** that extend the plate to the cervical spine and two **feet** that hold the rod-to-plate connection. The four sizes vary in length and width of the plate legs. The body and feet of the plate remain constant.



Midline Occiput Plate

Part #	Description
48551044	Midline Occiput Plate, Small
48551045	Midline Occiput Plate, Medium
48551046	Midline Occiput Plate, Large
48551047	Midline Occiput Plate, Large Long
48551048	Midline Occiput Plate, Mini



Midline Plate



Midline Plate, Mini

* Mini can be specially ordered

Implant Design Features

Plate Body

- Allows five points of fixation to the occiput, three points on the midline*, and two additional points at the superior lateral end
- Can be contoured in all planes to adapt to varying occipital anatomy
- Curved shape at the top of the plate is designed to be positioned just beneath the external occipital protuberance, and is constant for all sizes
- Hole configuration is the same for all sizes, allowing for interchangeability if needed.



Plate Legs

- Extend the plate to the cervical spine to help reduce rod contouring
- Length and width vary (See Table 1)
- Width of the legs should match the width of the cervical fixation points as best as possible to minimize rod contouring
- Can be contoured in all planes to adapt to varying occipital anatomy


Table 1

Part #	Description	Length	Width	# Holes
48551044	Midline Occiput Plate, Small	30	37	5
48551045	Midline Occiput Plate, Medium	38	40	5
48551046	Midline Occiput Plate, Large	49	45	5
48551047	Midline Occiput Plate, Large Long	49	49	5
48551048	Midline Occiput Plate, Mini*	27	33	4

*Special Order

Plate Feet

- Designed to sit above the ring of C1, in line with the cervical fixation points
- Pre-angled at 130° (cannot be contoured further)
- Polyaxial tulip head connects the rod to the plate, and offers 26° of variability in each direction



*The Mini plate offers two points of midline fixation

Midline Occiput Plate Placement

1. Select Plate Size

The proper Occiput Plates should be selected allowing application immediately below the level of theinion. Size selection should attempt to optimize the relationship between the polyaxial heads of the plate and the cervical fixation point. The pre-angled feet of the plate are designed to sit above the ring of C1.



Midline Plate



Midline Plate, Mini


**Plate positioned using
Plate Holder
48500300**



2. Plate Contouring

To contour the Midline Occipital Plate use the **Occiput Plate Bending Irons**. The **Occiput Plate Bending Irons** can be used in a variety of positions. Plates should be contoured to better accommodate the patient's occipital anatomy.

Warning: It is important to avoid excessive plate bending, especially over the midline screw holes.



The midline plate can be bent in both the coronal and sagittal planes.

To contour plate to the patient's occipital anatomy, first bend the **legs** of the plate (Image 1 and 2) in the sagittal plane.



Image 1



Image 2

To bend plate legs, place the **mouth of one bender** over the **superior portion of the plate leg** and the **legs of the bender** over the adjacent leg of the plate. Benders are marked “TO BEND LEGS ONLY” (Image 3).

Warning: Do not bend the angled foot of the plate as excessive bends can effect the integrity of the construct and/or damage the polyaxial tulip head.



Image 3

Continue to produce a symmetrical bend on the adjacent side, if necessary (Image 2).



Image 2

To create a **sagittal bend** of the plate body, place the **head of the bender** on the **head of the plate**. Place the **mouth** of the other bender on the **inferior portion of the plate body** (Image 4).



Image 4

To create a **coronal bend**, hold both benders with the **mouth** facing upwards. Benders should align (Image 5) to contour the plate down the midline. Bend of the plate body should align with final contour of the plate legs.



Image 5

Tip: Holding the far ends of the bending iron helps to contour the plate appropriately.

It is important to avoid point bending (Image 6) as excessive bending can damage the plate resulting in plate breakage.

Warning: Do not repeatedly bend and debend the plate. Care should be taken to not make extreme bends so as to avoid stress concentrations on the plate.



Image 6

Plate Placement

Use the **Plate Holder** to hold plate in position on the occiput while implanting the bone screws.



Plate Holder
48500300



Plate Holder on bone

Screw Placement

- a) Use the **Awl** to penetrate the cortex and mark the bone screw entry point.
- b) **Drill** depth should be determined based on approximating the thickness of the occiput from preoperative imaging and based on the anatomy. Drilling should be performed through the Occiput Plate using the appropriate Drill Guide. Two-sided Occiput Fixed Drill Guides are available in depths of 6 & 8mm, 10 & 12mm, 14 & 16mm and can be used with all diameter screws.




Awl
48560010

Occiput Fixed Drill Guide

Part #	Description
48561044	Occiput Fixed Drill Guide, 6 & 8mm
48561045	Occiput Fixed Drill Guide, 10 & 12mm
48561046	Occiput Fixed Drill Guide, 14 & 16mm



Note: It is recommended to Drill as close to full thickness into the occiput as can be safely performed. Drill depth should be determined based on approximating the thickness of the occiput from preoperative imaging and based on patient anatomy. Drilling should be performed through the Occiput Plate with the Drill Guide set. The Drill Guide should be long enough to accommodate the fact that the plate might not be flush to the bone.



The Occiput Drill Bits are available in three diameters and are color-coded to match the corresponding screw diameters:

- 2.5mm Drill for 3.5mm Screw
- 3.0mm Drill for 4.0mm Screw
- 3.5mm Drill for 4.5mm Screw

Screw Diameter	Corresponding Drill	Corresponding Tap
3.5mm 485543(06-16)* 	2.5mm Drill 48565323 	Occiput Set Tap, 3.5mm 48561053 
4.0mm 485544(06-16)* 	3.0mm Drill 48565423 	Occiput Set Tap, 4.0mm 48561054 
4.5mm 485545(06-16)* 	3.5mm Drill 48565523 	Occiput Set Tap, 4.5mm 48561055 

*2mm increments


Drill Bits can be attached to the **Short Quick Release Handle** for manual drilling. Alternatively, you can attach the Drill Bits to a power drill using either the **Flexible Shaft** or the **Straight Shaft**.



Warning: Be sure to use the corresponding Drill and Tap. Using a smaller diameter Drill can result in Tap breakage. Tapping in a hole too small can result in a sheared tap.


- c) Attach the appropriate sized **Occiput Set Tap** (color-coded to match the screw diameter) to the Short Quick Release Handle, and be sure to tap the entire length of the hole. Measure the length tapped to ensure appropriate screw length.

d) Insert the bone screw using the **Short Driver**, **Angled Driver** or the **Hexdriver**. No matter which driver is used to initially insert the bone screw, it is recommended that the Short Driver be used for final tightening. The maximum number of screws should be used as allowed by the anatomy.



3. Contour Rod

Contour the rod as required, using the Rod/Plate Bender. For a more acute angle, you may also bend the rods using the small holes of the Occiput Plate Bending Irons.



**Occiput Plate
Bending Irons**
48561049




Lay the rod into the tulip heads of the Occiput Plate and the cervical fixation points. The rod is fixed to the Occiput Plate with the standard blocker. Insert the blockers with the Hexdriver, utilizing the Persuader and/or Insertion Tube as necessary. Final tighten all blockers to 3Nm. The Rod must extend through the polyaxial head for proper fixation.


Use the Anti-Torque Key and Torque Wrench on the polyaxial heads. The standard blocker, which should be tightened by the Torque Wrench to 3Nm, is used to secure the plate-to-rod assembly.



Hexdriver inserts blocker into polyaxial tulip.



Final tighten blocker to 3Nm.



Blocker insertion with Hexdriver.

Bilateral Occiput Plate Placement

The OASYS® system also features bilateral plates that can be cut and contoured as required to the patient's anatomy. Both straight and pre-contoured Occiput Plates are available in 100° and 130° foot angles.

The Occiput Plates feature four round anchor holes and one slot, as well as a foot, which connects the plates into a rod construct in the cervical spine. The angle minimizes the amount of rod contouring necessary in the sagittal plane to fit the implants.




Bilateral Occiput Plates

Part #	Description
48551040	130° Occiput Plate
48551041L	130° Pre-contoured Occiput Plate, Left
48551041R	130° Pre-contoured Occiput Plate, Right
48551042	100° Occiput Plate
48551043L	100° Pre-contoured Occiput Plate, Left
48551043R	100° Pre-contoured Occiput Plate, Right

Both the 100° and 130° plates are offered in a “left” and “right” pre-contoured version. The plates have approximately a 30° radius, and cannot be contoured further in the frontal plane using the Rod/Plate Bender.

To prepare the Rod/Plate Cutter for multilevel plate insertion, move the handles into the “open” position as indicated with the laser-marked arrows. Insert the plate into the slot in the head of the Rod/Plate Cutter from the side that contains the laser marking (arrow), and position it so that the bending groove between holes lines up with the end of the small gauge block situated above the insertion slot. This position allows the plate to be cut through the area between holes. Compress the handles as indicated by the arrow to shorten the plate.



**Rod/Plate Cutter
48560018**

Note: If plate is overcontoured in lateral direction, replace with new implant.

Using the **Bilateral Occiput Plate Bend Template**, determine the desired position and contour of the implant. The side of the Rod/Plate Bender without the round bending iron is used to bend the plate in the frontal plane. Slide the plate in the undercut portion of the top bending arms and position the bending tip between the plate holes. Once bent in the frontal plane, the plate should never be reverted to its original shape as it would compromise its strength.



Rod/Plate Bender
48560019

The opposite side of the Rod/Plate Bender with the round bending iron is used to bend the plate in the sagittal plane. It is recommended that the plate be bent gradually from both sides and to avoid excessive contouring, which may compromise the plate material integrity. Bending the plate in the area between the screw holes minimizes deformation of the screw holes and allows proper screw fit.



Sagittal bending with Rod/Plate Bender

Note: The plate should always be bent in the frontal plane first. The plate should not be contoured more than 8° per segment. If plate is damaged or overbent, replace implant.

When the plate has been contoured to its desired shape and placed in position with the Plate Holder, use it as a guide to mark the entry points for the pilot holes with the Awl. As an alternative, the blue plate template can also be used for this purpose. Remove the plate or the plate template, and proceed with drilling the holes, using the Drill Guide and the appropriate size Drill (3.5, 4.0 or 4.5mm, depending on the diameter of the screw to be used). Drill depth should be determined based on approximating the thickness of the occiput from preoperative imaging and based on the anatomy. Drilling should be performed through the Occiput Plate with the Drill Guide set to the appropriate length. The hole should then be tapped.




Plate Holder and Awl

To connect the Occiput Plate to a cervical construct, it is recommended that the plate be pre-assembled onto the rod prior to implanting the plate. The standard blocker, which should be tightened by the Torque Wrench to 3Nm, is used to secure the plate-to-rod assembly. Connect the plate to the rod prior to putting in the bone screws.

Position the final plate implant, and insert and final tighten the occiput bone screws with the Angled Driver, Hexdriver or Short Driver.



Insert blocker into plate with Hexdriver



For increased stiffness and strength, an added Vitallium[®] rod option is available.

Part #	Description
48553080	3.5mm x 80mm Vitallium [®] Rod
48553120	3.5mm x 120mm Vitallium [®] Rod
48553240	3.5mm x 240mm Vitallium [®] Rod
48553350	3.5mm x 350mm Vitallium [®] Rod



Vitallium[®] Rod

Note: For increased construct strength, at least two bone screws should be placed in any of the four anchor holes. A bone screw should also be inserted into the slot at the base of the plate for a total of at least three screws in each plate. Crosslinks are recommended to increase construct stability and stiffness.

Implants

Part #	Description
--------	-------------

Biased Angle Screws

48552310	3.5mm x 10mm Biased Angle Screw
48552312	3.5mm x 12mm Biased Angle Screw
48552314	3.5mm x 14mm Biased Angle Screw
48552316	3.5mm x 16mm Biased Angle Screw
48552318	3.5mm x 18mm Biased Angle Screw
48552320	3.5mm x 20mm Biased Angle Screw
48552322	3.5mm x 22mm Biased Angle Screw
48552324	3.5mm x 24mm Biased Angle Screw
48552326	3.5mm x 26mm Biased Angle Screw
48552328	3.5mm x 28mm Biased Angle Screw
48552330	3.5mm x 30mm Biased Angle Screw
48552332	3.5mm x 32mm Biased Angle Screw
48552334	3.5mm x 34mm Biased Angle Screw
48552336	3.5mm x 36mm Biased Angle Screw
48552338	3.5mm x 38mm Biased Angle Screw
48552340	3.5mm x 40mm Biased Angle Screw
48552342	3.5mm x 42mm Biased Angle Screw
48552344	3.5mm x 44mm Biased Angle Screw
48552348	3.5mm x 48mm Biased Angle Screw
48552350	3.5mm x 50mm Biased Angle Screw
48552352	3.5mm x 52mm Biased Angle Screw
48552354	3.5mm x 54mm Biased Angle Screw

48552410	4.0mm x 10mm Biased Angle Screw
48552412	4.0mm x 12mm Biased Angle Screw
48552414	4.0mm x 14mm Biased Angle Screw
48552416	4.0mm x 16mm Biased Angle Screw
48552418	4.0mm x 18mm Biased Angle Screw
48552420	4.0mm x 20mm Biased Angle Screw
48552422	4.0mm x 22mm Biased Angle Screw
48552424	4.0mm x 24mm Biased Angle Screw
48552426	4.0mm x 26mm Biased Angle Screw
48552428	4.0mm x 28mm Biased Angle Screw
48552430	4.0mm x 30mm Biased Angle Screw
48552432	4.0mm x 32mm Biased Angle Screw
48552434	4.0mm x 34mm Biased Angle Screw
48552436	4.0mm x 36mm Biased Angle Screw
48552438	4.0mm x 38mm Biased Angle Screw
48552440	4.0mm x 40mm Biased Angle Screw
48552442	4.0mm x 42mm Biased Angle Screw
48552446	4.0mm x 46mm Biased Angle Screw
48552448	4.0mm x 48mm Biased Angle Screw
48552450	4.0mm x 50mm Biased Angle Screw

Part #	Description
--------	-------------

Medial Biased Angle Screws

48555320	3.5mm x 20mm Medial Biased Angle Screw
48555322	3.5mm x 22mm Medial Biased Angle Screw
48555324	3.5mm x 24mm Medial Biased Angle Screw
48555328	3.5mm x 28mm Medial Biased Angle Screw
48555332	3.5mm x 32mm Medial Biased Angle Screw
48555336	3.5mm x 36mm Medial Biased Angle Screw
48555340	3.5mm x 40mm Medial Biased Angle Screw

48555420	4.0mm x 20mm Medial Biased Angle Screw
48555422	4.0mm x 22mm Medial Biased Angle Screw
48555424	4.0mm x 24mm Medial Biased Angle Screw
48555428	4.0mm x 28mm Medial Biased Angle Screw
48555432	4.0mm x 32mm Medial Biased Angle Screw
48555436	4.0mm x 36mm Medial Biased Angle Screw
48555440	4.0mm x 40mm Medial Biased Angle Screw
48555444	4.0mm x 44mm Medial Biased Angle Screw
48555448	4.0mm x 48mm Medial Biased Angle Screw
48555452	4.0mm x 52mm Medial Biased Angle Screw

Smooth Shank Biased Angle Screws

48556322	3.5mm x 22mm Smooth Shank Biased Angle Screw
48556324	3.5mm x 24mm Smooth Shank Biased Angle Screw
48556326	3.5mm x 26mm Smooth Shank Biased Angle Screw
48556328	3.5mm x 28mm Smooth Shank Biased Angle Screw
48556330	3.5mm x 30mm Smooth Shank Biased Angle Screw
48556332	3.5mm x 32mm Smooth Shank Biased Angle Screw
48556334	3.5mm x 34mm Smooth Shank Biased Angle Screw
48556336	3.5mm x 36mm Smooth Shank Biased Angle Screw
48556338	3.5mm x 38mm Smooth Shank Biased Angle Screw
48556340	3.5mm x 40mm Smooth Shank Biased Angle Screw

Cancellous Biased Angle Screws

48557310	3.5mm x 10mm Cancellous Biased Angle Screw
48557312	3.5mm x 12mm Cancellous Biased Angle Screw
48557314	3.5mm x 14mm Cancellous Biased Angle Screw
48557316	3.5mm x 16mm Cancellous Biased Angle Screw
48557318	3.5mm x 18mm Cancellous Biased Angle Screw
48557320	3.5mm x 20mm Cancellous Biased Angle Screw
48557322	3.5mm x 22mm Cancellous Biased Angle Screw
48557324	3.5mm x 24mm Cancellous Biased Angle Screw

48557410	4.0mm x 10mm Cancellous Biased Angle Screw
48557412	4.0mm x 12mm Cancellous Biased Angle Screw
48557414	4.0mm x 14mm Cancellous Biased Angle Screw
48557416	4.0mm x 16mm Cancellous Biased Angle Screw
48557418	4.0mm x 18mm Cancellous Biased Angle Screw
48557420	4.0mm x 20mm Cancellous Biased Angle Screw
48557422	4.0mm x 22mm Cancellous Biased Angle Screw
48557424	4.0mm x 24mm Cancellous Biased Angle Screw

Part #	Description
--------	-------------

Non-Biased Angle Polyaxial Screws



48558310	3.5mm x 10mm Non-Biased Angle Screw
48558312	3.5mm x 12mm Non-Biased Angle Screw
48558314	3.5mm x 14mm Non-Biased Angle Screw
48558316	3.5mm x 16mm Non-Biased Angle Screw
48558318	3.5mm x 18mm Non-Biased Angle Screw
48558320	3.5mm x 20mm Non-Biased Angle Screw



48558410	4.0mm x 10mm Non-Biased Angle Screw
48558412	4.0mm x 12mm Non-Biased Angle Screw
48558414	4.0mm x 14mm Non-Biased Angle Screw
48558416	4.0mm x 16mm Non-Biased Angle Screw
48558418	4.0mm x 18mm Non-Biased Angle Screw
48558420	4.0mm x 20mm Non-Biased Angle Screw



48558520	4.5mm x 20mm Non-Biased Angle Screw
48558524	4.5mm x 24mm Non-Biased Angle Screw
48558528	4.5mm x 28mm Non-Biased Angle Screw
48558532	4.5mm x 32mm Non-Biased Angle Screw
48558536	4.5mm x 36mm Non-Biased Angle Screw
48558540	4.5mm x 40mm Non-Biased Angle Screw
48558544	4.5mm x 44mm Non-Biased Angle Screw
48558548	4.5mm x 48mm Non-Biased Angle Screw
48558552	4.5mm x 52mm Non-Biased Angle Screw

Rods

48553080	3.5mm x 80mm Vitallium [®] Rod
48553120	3.5mm x 120mm Vitallium [®] Rod
48553240	3.5mm x 240mm Vitallium [®] Rod
48553350	3.5mm x 350mm Vitallium [®] Rod

48552025	3.5mm x 25mm Titanium Alloy Rod
48552030	3.5mm x 30mm Titanium Alloy Rod
48552040	3.5mm x 40mm Titanium Alloy Rod
48552050	3.5mm x 50mm Titanium Alloy Rod
48552060	3.5mm x 60mm Titanium Alloy Rod
48552070	3.5mm x 70mm Titanium Alloy Rod
48552080	3.5mm x 80mm Titanium Alloy Rod
48552120	3.5mm x 120mm Titanium Alloy Rod
48552240	3.5mm x 240mm Titanium Alloy Rod
48551240	3.5mm x 240mm Titanium Rod, CP
48551350	3.5mm x 350mm Titanium Alloy Rod

Part #	Description
--------	-------------

Connectors



48551070	Transverse Connector, 80mm
48551071	Transverse Connector, 60mm
48551072	Transverse Connector, 40mm
48551073	Cross Connector Plate, 24mm
48551074	Cross Connector Plate, 32mm
48551075	Cross Connector Plate, 40mm
48551080	Offset Connector, 20mm
48551081	Offset Connector, 12mm
48551085	3.5mm to 3.5mm Axial Connector
48551084	3.5mm to 4.5mm Axial Connector
48551086	3.5mm to 5.5mm Axial Connector
48551087	3.5mm to 6.0mm Axial Connector
48551088	3.5mm to 3.5mm Parallel Rod-to-Rod Connector
48551091	3.5mm to 4.5mm Parallel Rod-to-Rod Connector
48551089	3.5mm to 5.5mm Parallel Rod-to-Rod Connector
48551090	3.5mm to 6.0mm Parallel Rod-to-Rod Connector



48551094	Saddle Connector
----------	------------------

48551000	Blocker
----------	---------

48551005	Set Screw
----------	-----------

Bone Screws



48554306	3.5mm x 6mm Bone Screw
48554308	3.5mm x 8mm Bone Screw
48554310	3.5mm x 10mm Bone Screw
48554312	3.5mm x 12mm Bone Screw
48554314	3.5mm x 14mm Bone Screw
48554316	3.5mm x 16mm Bone Screw



48554406	4.0mm x 6mm Bone Screw
48554408	4.0mm x 8mm Bone Screw
48554410	4.0mm x 10mm Bone Screw
48554412	4.0mm x 12mm Bone Screw
48554414	4.0mm x 14mm Bone Screw
48554416	4.0mm x 16mm Bone Screw














48554506	4.5 x 6mm Bone Screw
48554508	4.5 x 8mm Bone Screw
48554510	4.5 x 10mm Bone Screw
48554512	4.5 x 12mm Bone Screw
48554514	4.5 x 14mm Bone Screw
48554516	4.5 x 16mm Bone Screw

Part #	Description
--------	-------------

Hooks

	48551049 Hook, Standard, 3.5mm
	48551050 Hook, Standard Short, 5.0mm
	48551053 Hook, Standard, 6.5mm
	48551055 Hook, Standard Tall, 8.0mm
	48551060 Hook, Offset Right
	48551065 Hook, Offset Left

Occiput Fixation

	48551044 Midline Occiput Plate, Small
	48551045 Midline Occiput Plate, Medium
	48551046 Midline Occiput Plate, Large
	48551047 Midline Occiput Plate, Large Long
	48551048 Midline Occiput Plate, Mini
	48551040 130° Occiput Plate
	48551041L 130° Pre-contoured Occiput Plate, Left
	48551041R 130° Pre-contoured Occiput Plate, Right
	48551042 100° Occiput Plate
	48551043L 100° Pre-contoured Occiput Plate, Left
	48551043R 100° Pre-contoured Occiput Plate, Right

Instruments

Part #	Description
--------	-------------

Trays



48561001	Container 1, Polyaxial Screws
48561001A	Cancellous Screw Caddy
48561001B	Non-Biased Angle Screw Caddy
48561002	Container 2, Implants & Screw Prep
48561003	Container 3, General Instruments
48560002	Rod Cutter Tray
48561004	Container 4, Occiput Fixation Tray

Instruments



48560010	Awl
48561011	Adjustable Drill Guide
48561034	Fixed Drill Guide, 12mm
48561035	Fixed Drill Guide, 14mm
48560012	Quick Release Handle
48561052	Short Quick Release Handle
48561051	Flexible Shaft
48561056	Straight Shaft

Part #	Description
--------	-------------



48560323	2.5mm Drill for 3.5mm Screw
48560423	3.0mm Drill for 4.0 & 4.5mm Screws
48565323	2.5mm Drill for 3.5mm Bone Screw
48565423	3.0mm Drill for 4.0mm Bone Screw
48565523	3.5mm Drill for 4.5mm Bone Screw
48561044	Occiput Fixed Drill Guide, 6 & 8mm
48561045	Occiput Fixed Drill Guide, 10 & 12mm
48561046	Occiput Fixed Drill Guide, 14 & 16mm
48560314	3.5mm Tap
48560414	4.0mm Tap
48560514	4.5mm Tap
48561916	3.5mm Cancellous Tap
48561917	4.0mm Cancellous Tap
48561914	Tap Sleeve, 3.5 & 4.0mm
48561915	Tap Sleeve, 4.5mm
48561053	Occiput Set Tap, 3.5mm
48561054	Occiput Set Tap, 4.0mm

Part #	Description	Part #	Description
 48561055	Occiput Set Tap, 4.5mm	 48560022	<i>In-situ</i> Bender, Right
 48560015	Depth Gauge	 48560023	<i>In-situ</i> Bender, Left
 48561036	Ball Tip Probe	 48561049	Occiput Plate Bending Irons (2)
 48561115	Pedicule Probe	 48561024	Locking Persuader
 48562016	Polyaxial Screwdriver	 48560024	Persuader
 48561037	Poly Adjustment Driver	 48561042	Hexdriver
 48561043	Angled Driver	 48561026	Compressor
 48561050	Short Driver	 48561041	Distractor
 48560017	Rod Template	 48560027	Insertion Tube/Anti-torque Key
 48560050	Bilateral Occiput Plate Bend Template	 48560028	Torque Wrench
 48560018	Rod/Plate Cutter	 48560029	Connector Clip Inserter
 48561038	Small Rod Cutter	 48560129	Connector Clip Remover
 48560019	Rod/Plate Bender	 48561071	Connector Plate Bender
 48560020	Rod Forceps	 48561073	Locking Nut Socket, Connector Plate
 48561039	Rod Rotation Forceps	 48500300	Plate Holder
 48560021	Screw Head Adjuster	 48560032	Hook Preparer
		 48560033	Hook Forceps

STRYKER SPINE OASYS[®] SYSTEM

DESCRIPTION

The STRYKER Spine OASYS[®] System is intended for use as an aid in spine fusion. It consists of screws, hooks, plates, rods and connectors. These components are available in a variety of lengths in order to accommodate individual patient physiology and pathology and to facilitate posterior stabilization of the spine.

INDICATIONS

When intended to promote fusion of the cervical spine and occipito-cervico-thoracic junction (Occiput-T3), the STRYKER Spine OASYS[®] System is intended for:

- Degenerative Disc Disease (as defined by neck and back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies)
- Spondylolisthesis
- Spinal Stenosis
- Fracture/Dislocation
- Atlanto/axial fracture with instability
- Occipitocervical dislocation
- Revision of previous cervical spine surgery
- Tumors

When used with the occipital plate, the bone screws are limited to occipital fixation only. The bone screws are not intended to be used in the cervical spine.

The use of the polyaxial screws is limited to placement in the upper thoracic spine (T1-T3). They are not intended to be placed in the cervical spine.

The hooks and rods are also intended to provide stabilization to promote fusion following reduction of fracture/dislocation or trauma in the cervical/upper thoracic (C1-T3) spine.

The Stryker Spine OASYS[®] System can also be linked to the Xia[®] System, SR90D System and Xia[®] 4.5 System via the rod-to-rod connectors and polyaxial screws of Xia[®] II and Xia[®] 3 Systems via the saddle connector.

CONTRAINDICATIONS

Contraindications may be relative or absolute. The choice of a particular device must be carefully weighed against the patient's overall evaluation. Circumstances listed below may reduce the chances of a successful outcome:

- Any abnormality present, which affects the normal process of bone remodeling including, but not limited to, severe osteoporosis involving the spine, bone

absorption, osteopenia, active infection at the site or certain metabolic disorders affecting osteogenesis.

- Insufficient quality or quantity of bone, which would inhibit rigid device fixation.
- Previous history of infection.
- Excessive local inflammation.
- Open wounds.
- Any neuromuscular deficit, which places an unusually heavy load on the device during the healing period.
- Obesity. An overweight or obese patient can produce loads on the spinal system, which can lead to failure of the fixation of the device or to failure of the device itself.
- Patients having inadequate tissue coverage of the operative site.
- Pregnancy.
- A condition of senility, mental illness, or substance abuse. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the implant, leading to failure or other complications.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests must be made prior to material selection or implantation.
- Other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of, congenital abnormalities, elevation of sedimentation rate unexplained by other diseases, elevation of white blood cell count (WBC), or marked left shift in the WBC differential count.

These contraindications can be relative or absolute and must be taken into account by the physician when making his decision. The above list is not exhaustive.

INFORMATION FOR THE PATIENT

The surgeon must discuss all physical and psychological limitations inherent to the use of the device with the patient. This includes the rehabilitation regimen, physical therapy, and wearing an appropriate orthosis as prescribed by the physician. Particular discussion must be directed to the issues of premature weightbearing, activity levels, and the necessity for periodic medical follow-up.

The surgeon must warn the patient of the surgical risks and the patient must be made aware of possible adverse effects. The surgeon must warn the patient that the device cannot and does not replicate the flexibility, strength, reliability or durability of normal healthy bone, that the implant can

break or become damaged as a result of strenuous activity or trauma, and that the device may need to be replaced in the future. If the patient is involved in an occupation or activity which applies inordinate stress upon the implant (e.g., substantial walking, running, lifting, or muscle strain) the surgeon must advise the patient that resultant forces can cause failure of the device. Patients who smoke have been shown to have an increased incidence of non-unions. Surgeons must advise patients of this fact and warn of the potential consequences. For diseased patients with degenerative disease, the progression of degenerative disease may be so advanced at the time of implantation that it may substantially decrease the expected useful life of the appliance. In such cases, orthopaedic devices may be considered only as a delaying technique or to provide temporary relief.

ADVERSE EFFECTS

- While the expected life of spinal implant components is difficult to estimate, it is finite. These components are made of foreign materials, which are placed within the body for the potential fusion of the spine and reduction of pain. However, due to the many biological, mechanical and physicochemical factors, which affect these devices but cannot be evaluated in vivo, the components cannot be expected to indefinitely withstand the activity level and loads of normal healthy bone.
- Bending, disassembly or fracture of any or all implant components.
- Fatigue fracture of spinal fixation devices, including screws, rods, and hooks has occurred.
- Pain, discomfort, or abnormal sensations due to the presence of the device.
- Pressure on skin from components where inadequate tissue coverage exists over the implant, with the potential extrusion through the skin.
- Dural leak requiring surgical repair. This risk is related to the surgical procedure. The intended use of the device does not require it to be close to the dura.
- Cessation of growth of the fused portion of the spine.
- Loss of proper spinal curvature, correction, height and/or reduction.
- Delayed Union or Nonunion: Internal fixation appliances are load sharing devices which are used to obtain alignment until normal healing occurs. In the event that healing is delayed, does not occur, or failure to immobilize the delayed/nonunion results, the implant will be subject to excessive and

repeated stresses which can eventually cause loosening, bending or fatigue fracture. The degree or success of union, loads produced by weight bearing, and activity levels will, among other conditions, dictate the longevity of the implant. If a nonunion develops or if the implants loosen, bend or break, the device(s) must be revised or removed immediately before serious injury occurs.

- Loosening of spinal fixation implants can occur. Early mechanical loosening may result from inadequate initial fixation, latent infection, premature loading of the prosthesis or trauma. Late loosening may result from trauma, infection, biological complications or mechanical problems, with the subsequent possibility of bone erosion, migration and/or pain.
- Cervical spine procedures may be associated with vascular and neural complications such as arterial injury or mechanical compromise, cord contusion and damage, peripheral nerve compromise and damage, including but not limited to peripheral paralysis, sensory disorders, vascular disorders, loss or disturbance of bladder and bowel functions.
- Serious complications may be associated with any surgery. These complications include, but are not limited to: genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; bursitis, hemorrhage, myocardial infarction, infection, paralysis or death.
- Neurological, vascular, or soft tissue damage due directly to the unstable nature of the fracture, or to surgical trauma.
- Inappropriate or improper surgical placement of this device may cause distraction or stress shielding of the graft or fusion mass. This may contribute to failure of an adequate fusion mass to form.
- Decrease in bone density due to stress shielding.
- Intraoperative fissure, fracture, or perforation of the spine can occur due to implantation of the components. Postoperative fracture of bone graft, the intervertebral body, pedicle, and/or lateral mass above and/or below the level of surgery can occur due to trauma, the presence of defects, or poor bone stock.
- Adverse effects may necessitate reoperation or revision.
- The surgeon must warn the patient of these adverse effects as deemed necessary.

REMOVAL OF IMPLANTS

These implants are temporary internal fixation devices designed to stabilize the operative site during the normal healing process. After healing occurs, these devices serve no functional purpose and can be removed. Removal may also be recommended in other cases, such as:

- Corrosion with a painful reaction
- Migration of the implant, with subsequent pain and/or neurological, articular or soft tissue lesions
- Pain or abnormal sensations due to the presence of the implants
- Infection or inflammatory reactions
- Reduction in bone density due to the different distribution of mechanical and physiological stresses and strains
- Bone growth restraint due to the presence of the implants
- Failure or mobilization of the implant

Instruments are provided by STRYKER Spine to be used to remove the implants. Any decision by a physician to remove the internal fixation device must take into consideration such factors as the risk to the patient of the additional surgical procedure as well as the difficulty of removal. Removal of an unloosened spinal screw may require the use of special instruments to disrupt the interface at the implant surface. This technique may require practice in the laboratory before being attempted clinically. Implant removal must be followed by adequate postoperative management to avoid fracture or re-fracture. Removal of the implant after fracture healing is recommended. Metallic implants can loosen, bend, fracture, corrode, migrate, cause pain or stress shield bone.

PRE-OPERATIVE PRECAUTIONS

Surgical Technique brochures may be requested from a distributor or from STRYKER Spine directly. Those using brochures published more than two years before the surgical intervention are advised to obtain an updated version.

STRYKER Spine devices can only be used by doctors who are fully familiar with the surgical technique required and who have been trained to this end. The doctor operating must take care not to use the instruments to exert inappropriate stress on the spine or the implants and must scrupulously comply with any operating procedure described in the surgical

technique provided by STRYKER Spine. For example, the forces exerted when repositioning an instrument in-situ must not be excessive as this is likely to cause injury to the patient.

To reduce the risks of breakage, care must be taken not to distort the implants or nick, hit or score them with the instruments unless otherwise specified by the applicable STRYKER Spine Surgical technique.

Extreme care must be taken when the instruments are used near vital organs, nerves or vessels.

Unless otherwise specified on the label, the instruments may be reused after decontamination, cleaning and sterilization.

WARNINGS

Federal law restricts this device to sale by or on the order of a licensed physician. The safety and effectiveness of pedicle screw spinal systems have been established only for spinal conditions with significant mechanical instability requiring fusion with instrumentation. These conditions are significant mechanical instability of the thoracic, lumbar, and sacral spine secondary to degenerative spondylolisthesis with objective evidence of neurological impairment, fracture, dislocation, spinal tumor, and failed previous fusion (pseudoarthrosis). The safety and effectiveness of these devices for any other conditions are unknown.

PRECAUTIONS

The implantation of pedicle screw spinal systems must be performed only by experienced spinal surgeons with specific training in the use of this pedicle screw spinal system because this is a technically demanding procedure presenting a risk of serious injury to the patient. Based on the fatigue testing results, the physician/surgeon must consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact on the performance of the system.

The STRYKER Spine OASYS® System has not been evaluated for safety and compatibility in the MR environment. The STRYKER Spine OASYS® System has not been tested for heating or migration in the MR environment.

